

# Coherent Resonant Soft X-Ray Magnetic Scattering

Truth from speckle  
How magnets forget

Larry Sorensen  
Department of Physics  
University of Washington

How can we determine  
precisely what the ensemble  
of nanoscopic magnetic  
domains is doing during  
static and dynamic magnetic  
hysteresis in the presence of  
disorder?

# MY COWORKERS:

Michael Pierce	PhD	University of Washington
Paul Unwin	MS	
Elaine Chan	MS	
Robert Moore		
Bo Hu		
Phillip Geissbuhler		

Steve Kevan	University of Oregon
Josh Turner	

Jeffrey Kortright	Lawrence Berkeley Laboratory
Karine Chesnel	

Eric Fullerton	IBM Almaden
Olav Hellwig	

All of these experiments were done  
using beamline 9.0.1 at the ALS.  
Raw undulator x-rays!

# What is the Ideal Source for Coherent Soft X-Ray Experiments?

## Lessons from Lasers:

- 1) CW (continuous wave)
- 2) Modulation (knobs!)
- 3) Below damage threshold


(1) Quasi-continuous beam is required to sample the thermal fluctuations efficiently. The proposed Stanford pulsed diffraction-limited undulator will be very bad for dynamic light scattering studies. No one has ever done dynamic light scattering with a pulsed laser.

Repetition rate 500 MHz to 500 GHz

(2) Variable bandwidth---from  $10^{-2}$  to  $10^{-5}$ ---is required to provide the optimum longitudinal coherence length for different dynamic light scattering experiments.

The narrow bandwidth beams from such a machine will have complete transverse coherence, and the longitudinal coherence necessary to do dynamic light scattering studies without damaging the samples---again, this is in sharp contrast to the proposed Stanford source.

# COMPARE WITH LASER SPECKLE



## LASER SPECKLE

### Coherent flux in

Wavelength

Transverse correlation length

Longitudinal correlation length

$10^{16}$  photons/second

5320 Angstroms

2 mm (4,000 )

10 cm (200,000 )

## UNDULATOR SPECKLE

### Coherent flux in

Wavelength

Transverse correlation length

Longitudinal correlation length

$10^{12}$  photons/second

16 Angstroms

30 microns (19,000 )

640 Angstroms (40 )

### Scattered flux out

Transverse correlation length

Longitudinal correlation length

Average size of magnetic domains

Number of magnetic domains

$10^6$  to  $10^7$  photons/second

30 microns (19,000 )

25,600 Angstroms (1600 )

2,000 Angstroms

40,000

# DIFFERENT APPELLATIONS FOR DYNAMIC LIGHT SCATTERING

CXRD: coherent x-ray diffraction (CXD)

DLS: dynamic light scattering (SXDLS, XDLS)

FPI: Fabry-Perot Interferometry

IFS: intensity fluctuation spectroscopy (XIFS)

II: intensity interferometry

HBTI: Hanbury-Brown Twiss interferometry

HES: heterodyne spectroscopy

HOS: homodyne spectroscopy

LBS: light beating spectroscopy

OBS: optical beating spectroscopy

PCS: photon correlation spectroscopy (XPCS)

PDII: post detection intensity interferometry

QELS: quasielastic light scattering

# WORLD SPEED RECORDS:

Liquid crystal samples

1999 106 times faster

fast scintillator + PM

1 microsecond limit

3 microsecond times


2003 20 times faster

avalanche photodiode

50 nanosecond limit

300 nanosecond times






Speckle Patterns versus  
Applied Magnetic Field



# Photons at the ALS




raw flux =  $2 \times 10^{14}$  photons/sec

coherent flux =  $2 \times 10^{12}$  photons/sec

scattered flux =  $2 \times 10^7$  photons/sec



## 0.01% Bandwidth



# COHERENT FLUX COMPARISON:

## Proposed SXU

50 x ALS-2

damage threshold / 50

## Theoretical SXU

770 x ALS-2

damage threshold / 3.3

## Ultimate Ring

2600 x ALS-2

500 x PSXU NSLS-2

33 x TSXU NSLS-2

# RESOLUTION LIMITS:

Proposed SXU	1 microsec
Theoretical SXU	30 nanosec
Ultimate Ring	0.5 nanosec

# RECOMMENDATIONS:

## 1) Source Characteristics:

Quasi-CW and just below damage threshold

## 2) Photon Ecology:

Use all photon degrees of freedom:

- a) Momentum
- b) Energy
- c) Polarization
- d) Coherence

=> Fast modulation!

OPTICS:


Perfect Optics


No Optics

Just Right Optics



# Scan Undulator Gap





# DETECTORS:

Perfect CCD

Fast Parallel

Associated Correlators

COHERENCE:

Perfect Coherence  
versus


Partial Coherence

Want Variable Transverse  
and Longitudinal Coherence!

# Magnetic Speckle

Coherent illumination of many  
(randomly located) magnetic  
domains => speckle pattern

Moving even a single magnetic domain  
produces large changes in the speckle  
pattern.



By measuring these changes versus  
the applied magnetic field (history), we can  
determine the changes in the sample---both  
in time and space.



# MAGNETIC SPECKLE

Coherent x-rays => coherent illumination

Add up the scattering amplitudes from each magnetic domain.

Randomly positioned magnetic domains => speckle pattern

Moving even a single magnetic domain can produce a **large change** in the speckle pattern.

By measuring the changes in the speckle pattern, we can determine the changes in the sample---both in space and time.

# COHERENT SOFT X-RAY MAGNETIC SCATTERING

The magnetic domains in our sample act like many superimposed diffraction gratings. The photon scattering from these gratings produces our magnetic speckle patterns.

By measuring these magnetic speckle patterns, we can precisely determine very small changes in the magnetic domain arrangement.

Because our photon degeneracy parameter  $\delta$  is very small, the interference takes place one photon at a time---our  $\delta$  is only about  $10^{-3}$ .

Hanbury-Brown and Twiss first developed post detection intensity interferometry for optical and RF stellar interferometry.